

# FACE MASK & FILTERS STERILIZER



### Inner Sterilization Space:

198\*105\*30 mm

### Product size:

231\*127\*47 mm

### Weight:

250 g

### Connection:

USB Cable/micro spin Type/5.0Vdc, 1.0A

### Standards

CE FC RoHS

### Description

The GVS Face Mask Sterilizer is an efficient solution for sterilizing F31000 FFP3 masks, as well as the filters from the Elipse P3 range. A LED lamp inside the sterilizer uses UVC light, to kill viruses in minutes, without damaging the mask or filter. The sterilizer can also be used to give a mask a pleasant scent via the Aromatherapy option.

### Features

- ◆ High efficiency: UVC lamp rays with a wavelength of 253.7nm sterilize a mask/filter removing 99.9% of contaminants in just ten minutes per side.
- ◆ The aromatherapy option allows the user to add a pleasant scent to their mask.
- ◆ A compact, portable and stylish product.
- ◆ The sterilizer doesn't use heat, liquids or chemicals to sanitize masks - just UV Light.
- ◆ Fully CE approved and certified.
- ◆ The product was tested with GVS respirators and Elipse P3 filters.
- ◆ Can also be used for small objects.

| Product Code | Description          | Quantity   |
|--------------|----------------------|------------|
| SPM900       | Face Mask Sterilizer | 20 per box |

## Applications



Travel



Medical



At Home



At Work